


DH-II Rotary Blender



DNA blender and vertical blender are widely used in sample mixture of medicine, agronomy and bioengineering, substitute rock by hand, convenient to use and simple to operate, can be put into all kinds of centrifugal pipes at the same time once, adjustable rotate speed of 10-60 turns/min., continuous work in 24 hours in succession, which can be put into the refrigerator

HS-3 Upright Blender



This product is used in chemistry, immunology, biological mixing reaction extensively, especially suitable for the mixing reaction of immune magnetic pearl. It can be put into constant temperature or refrigerator. Easy to use, simple to operate, can be put all kinds of centrifuge tube at odds with adjustable rotate speed, allow to work in succession in 24 hours.

Main parameters

Rating voltage: 220V 50HZ

Motor: 12VDC 60r/min

Rotate speed: 10~60 r/min

Specification of clamping centrifuge tube: 1.5-30ml