


ELISA Plates



Cat. No.	Size	Description	Binding Capacity	Quantity
FEP10096	96 well	Fixed, Flat	High binding capacity, 300~400 ng IgG/cm ² , suitable as solid phase of protein with medium /large molecular weight (> 10kD)	10/pkg, 200/case
FEP11096	96 well	Detachable, Flat		10/pkg, 200/case
FSP11048	48 well	Detachable, Flat		10/pkg, 200/case
FSP11012	12 well	Strip		40/pkg, 1600/case
FSP11008	8 well	Strip		60/pkg, 2400/case
FEP20096	96 well	Fixed, Flat	Medium binding, 200~300 ng IgG/cm ² , suitable as solid phase of protein with large molecular weight (> 20kD)	10/pkg, 200/case
FEP21096	96 well	Detachable, Flat		10/pkg, 200/case
FSP21048	48 well	Detachable, Flat		10/pkg, 200/case
FSP21012	12 well	Strip		40/pkg, 1600/case
FSP21008	8 well	Strip		60/pkg, 2400/case

FEP30096	96 well	Fixed, Flat	Aminate , binding capacity of 100 ng IgG /cm ² . Suitable as solid phase of protein with small molecules	10/pkg, 200/case
FEP31096	96 well	Detachable, Flat		10/pkg, 200/case
FSP31048	48 well	Detachable, Flat		10/pkg, 200/case
FSP31012	12 well	Strip		40/pkg, 1600/case
FSP31008	8 well	Strip		60/pkg, 2400/case

ELISA Plate Frames

Cat. No.	Size	Description	Quantity
EPF41096	96 well Frame	12w x 8或8w x 12	10/pkg, 200/case
EPF41048	48 well Frame	12w x 4	10/pkg, 200/case