


# Mini Cordless Grinder & Pestles

## Mini Cordless Grinder

### Features & Benefits

- Cordless, Portable, Easy to Operate.
- Battery Operation (4 pcs of AA Batteries)
- Applicable for Usage with 1.5cc & 0.5cc Disposable Pestle

**Cat# CG**



## 1.5cc & 0.5cc Disposable Pestle

### Features & Benefits

- Two Options: 1.5cc & 0.5cc Disposable Pestle
- P.P. material, Autoclavable.
- RNase-free, DNase-free & Pyrogen-free
- Applicable for usage with Mini Cordless Grinder (#CG-4A) or Manual Separately
- Economic Package: 100's/bag, 10bags/case

#### Cat No.

**3000-1005:** 1.5cc Disposable Pestle, PP, 100's/bag, 10bags/case (Suggest to use with "SSI" #HT15NG microcentrifuge tube for better effect.)

**3000-1004:** cc Disposable Pestle, PP, 100's/bag, 10bags/case (Suggest to use with "SSI" #HT05NG microcentrifuge tube for better effect.)