

# XHF-D High-Speed Homogenizer



This two kinds of instruments carry out work by electron timing control and operation mini high speed motor with special stainless steel inside and outside club crushing knife. By tail-wagging of cutting knife, it goes on mix, emulsification, and crash to various materials.

Technical performance and index of XHF-D high speed dispersing device

Voltage: 220V

Power: 360W

Timing range: 5000-24000turns/min.

Timing range:1-99min. (computer control, setting at will)

Ambient temperature: -10<sup>0</sup>C-40<sup>0</sup>C

Weight of whole machine: 6.0Kg

Overall dimension: 320×220×430

Technical performance and index of XHF-I high-speed dispersing device

Voltage: 220V

Power: 360W

Timing range: 2800-28000turns/min.

Timing range:5 sec, 10 sec, 20 sec... 9-grade continuous timing in total

Ambient temperature: -10<sup>0</sup>C -40<sup>0</sup> C

Weight of whole machine: 6.0Kg

Overall dimension: 320×220×430

# DY89-II Electric Glass Homogenizer



This machine utilizes mechanism of grinding, carry on soft melting evenly treatment in moving plant tissues, and get the cytoplasm, mitochondria of the experiment need. Suitable for fields such as biology, medical science, agronomy, etc., have low-speed big moment and noiseless characteristics etc.

Speed	50-1500rpm Constantly adjustable
Wattage	120W
Dimension	10ml.20ml with one
Special orders	1ml,3ml.5ml.50ml

**DY89-I without rotary speed display , Others are same as parameters of DY89-II.**