

Electrophoresis Power Supply



DYY-6C

Model	DYY-6C
Volt	600
Voltage range	5 ~ 600V
Current	4 ~ 400mA
Power	240W
Programs	---
Output terminals	4



DYY-12C


Model	DYY-12C
Volt	5000
Voltage range	10 ~ 5000V
Current	2 ~ 200mA
Power	2~200W
Programs	9
Output terminals	2



DYY-10C


Model	DYY-10C
Volt	3000
Voltage range	10 ~ 3000V
Current	2 ~ 200mA
Power	5~200W
Programs	9
Output terminals	2

Electrophoresis Power Supply



DYY-7C

Model	DYY-7C
Volt	200
Voltage range	2 ~ 200V
Current	10 ~ 2000mA
Power	200W
Programs	---
Output terminals	4



DYY-8C

Model	DYY-8C
Volt	600
Voltage range	5 ~ 600V
Current	2 ~ 200mA
Power	120W
Programs	---
Output terminals	2