

# Balances

Patented Products from GoldTop Import & Export Co. Ltd



Model (WT)	1002	2002	3002	6002	10001	20001	30001	60001
Maximum Weight	100g	200g	300g	500g	1000g	2000g	3000g	5000g
Minimum Reading	0.01g				0.1g			
Tray Size	Diameter 30mm				Diameter 60mm			
Description	g, oz, ct, GN, dwt conversion; simple counting; percentage weighting; RS232; printing; weight value correction; AC/DC; linear adjustment; bottom hook weighting; built-in calibrating weights optional for real-time automatic calibration.				g, oz, ct, GN, dwt conversion; simple counting; percentage weighting; RS232; printing; weight value correction; AC/DC; linear adjustment; bottom hook weighting; built-in calibrating weights optional for real-time automatic calibration.			

Patented technology from China (Patent Application number 03356372.1) from GoldTop Import & Export Co. Ltd ([www.machineweigh.com](http://www.machineweigh.com))



Model	LA110	LA160	LA200	LB120	LB200	LB300
Capacity	110g	160g	200g	120g	200g	300g
Reliability	0.1 mg			1 mg		
Linear Error	<+ 0.2 mg			<+2 mg		
Reproducibility	<- 0.1 mg			<-1 mg		
Calibration way	External calibration with a standard weight					
Data Output	RS232 interface					
Tare Range	0-110g	0-160g	0-200g	0-120g	0-200g	0-300g
Pan Size	Ø785mm					
Dimensions	315x225x330					
Power Supply	220V AC					