

Basal Cell Culture Media

Product name	Catalogue number	500ml <i>Liquid</i>	10×1L <i>Powder</i>	1×10L <i>Powder</i>
BME With Earle's salts and L-glutamine. Without sodium bicarbonate. Used for mammalian diploid cell culture	MD102	Enquire price	Enquire price	Enquire price
DMEM (High glucose) With D-glucose(4500mg/L) and L-glutamine. Without sodium bicarbonate	MD201	Enquire price	Enquire price	Enquire price
DMEM (High glucose) With D-glucose (4500mg/L), L-glutamine ,sodium pyruvate (110mg/L). Without sodium bicarbonate	MD202	Enquire price	Enquire price	Enquire price
DMEM/F-12 With L-glutamine, pyridoxine hydrochloride and pyridoxal hydrochloride. Without HEPES and sodium bicarbonate	MD207	Enquire price	Enquire price	Enquire price
E199 With Earle's salts and L-glutamine. Without sodium bicarbonate	MD500	Enquire price	Enquire price	Enquire price
MEM With Earle's salts and L-glutamine. Without sodium bicarbonate	MD600	Enquire price	Enquire price	Enquire price
MEM (Can be autoclaved) With Earle's salts. Without L-glutamine and sodium bicarbonate	MD605	Enquire price	Enquire price	Enquire price
RPMI1640 With L-glutamine. Without sodium bicarbonate	MD800	Enquire price	Enquire price	Enquire price
Fischer's With L-glutamine. Without sodium bicarbonate	MD300	Enquire price	Enquire price	Enquire price