

Low-Serum Cell Culture Media

Product Name	Catalogue Numer	500ml Liquid	10×1L Powder	1×10L Powder
SXQ-1 The Vero cell use the low-serum cell culture medium by mini-carrier bioreactor, the serum dosage is no more than 3%	SLM100	Enquire price	Enquire price	Enquire price
199 SLM Completely animal component-free, substitute for the hydrolytic lactoprotein Earle's balance salt solution. Used for Chinese hamster kidney cell and chick embryo cell culturing	SLM110	Enquire price	Enquire price	Enquire price
199 SLM-BR The modified low-serum of 199 cell culture medium. Used for the low serum culture of the Vero cell by mini-carrier bioreactor	SLM111	Enquire price	Enquire price	Enquire price
MEM SLM The modified low-serum of MEM. Can used for the low serum culture of the BHK21 cell	SLM120	Enquire price	Enquire price	Enquire price