

## Magnetic Stirrers

The function of a stirrer is to agitate liquids for speeding up reactions or improving mixtures. Whether with heating or without heating, our magnetic stirrers are built for durability and dependability. Likewise the shafted stirrers with propeller are available to fit the most demanding applications including very viscous liquid stirring. All of our stirring motors are designed for hour after hour of continuous and suitable for unsupervised operation (a Xiangquan Chromtech product).



Model	Speed Range rpm	Max samples	Heating Range	Temperature control	Temperature Resolution	Disc Diameter
AM-3250A	100-1200	10L	Room temperature~300°C			150mm
AM-3250B	100-1200	10L	Room temperature~199°C	Digital	±2 °C	150mm
AM-3250C	100-1200	10L	Room temperature~199 °C		±2 °C	150mm
AM-5250A	100-1500	10L	Room temperature~380 °C			150mm
AM-5250B	100-1500	10L	Room temperature~380 °C	Digital	±2 °C	150mm
AM-5250C	100-1500	10L	Room temperature~380 °C	Digital	±2 °C	150mm