

Microcentrifuges

A product from Bertec Enterprise Co., Ltd



This **Micro Centrifuge** is **small**, **easy** to load and operate, **simple** to close the lid to operate and open it to stop, **durable** for strong hinge supports and latch lock, and practical for **complete parts** enclosed. Our unique design is **ULTRA QUITE**, it does not make any noise when in operation.

Only exchange your plugs, you can safely used in the voltage from 100V to 240V without worrying about the damage in wrong voltage.

There are three colors for your choice, black, green and ivory!!

Ultra Quite
CE Certified

Advantages

● **All Range Voltage**
Suitable: 100V-240V



● **Durable**

Never broken for latch lock and the hinge of the transparent cover.



● **Complete Parts**

With 6-tubes in 1.5cc/2cc rotor, 6 adaptors for 0.5cc tube, and 8-strip for 0.2cc tubes.



● **Safe**

Turn on/off, without opening the lid.



● **Small**

Only 6"x 6" and 828g



● **Easy Change Rotor**

Very easy to change the rotors within 30 seconds



PC-100 Specifications

Speed	6,200 rpm (+20%)
Force(g)	2,000 x g
	6x1.5 to 2.2 ml tubes
Maximum Volume	2 Strips (8x 200ul tubes)
Power Required	100~ 240V, 260mA, 50/60 Hz
Dimensions	140Wx125Dx115H (mm)