


pH Meters

The meters for pH, ion concentrations and conductivity meet all the requirements of modern field and laboratory measurements (a Xiangquan Chromtech product).



Model	Measuring ranges			Resolution		
	pH	Potential	Temp.	pH	Potential	Temp.
PHS-25J	0.00~14.00	-1400~1400		0.01	1	
PHS-2F	0.00~14.00	-1400~1400		0.01	1	
PHB-4 mobile	0.00~14.00	-1400~1400		0.01	1	
PHBJ-260 mobile	0.00~14.00	-1800~1800	0.0~60.0	0.01	0.1	0.1
Model	Measuring ranges			Resolution		
	pH	Potential	Temperature	pH	Potential	Temperature
PHREx-2 mobile	0~14.0			0.1		
PHS-3B	0.00~14.00	-1999~1999	0.0~100.0	0.01	1	0.1
PHS-3C	0.00~14.00	-1999~1999		0.01	1	
PHSJ-3F	-2.00~18.00	-1999~1999	0.0~100.0	0.01	1	0.1
PHS-3E	-2.00~18.00	-1999~1999	0~99.9	0.01	1	0.1
PHS-3D	0.00~14.00	-1999~1999		0.01	1	
PHSJ-4A Lab	0.00~14.00	-1999~1999	-5.0~105.0	0.001	0.1	0.1
PHSJ-5 Lab	-2.00~18.00	-1999~1999	-5.0~105.0	0.001	0.1	0.1