

Shakers

Shakers are designed for applications that require precise control with speed and mode of operation. Small, quiet shaker ensures long life with ideal swivel motion. A digital display makes it possible to read the speed, timer function and operating modes (a Xiangquan Chromtech product).



Model	Speed Range(rpm)	Variation Range(mm)	Capacity(kg)	Time adjustable (min)	Temperature
ARCS-50	0~350	10	5	3~60	
ARS-50	0~150	3	4	3~60	
AVS-50	0~2500	3	1	3~60	
AVS-100	0~2500	4.5	1	3~60	
HMS-3000	0~3000	3	4	0~60	
HMS-310	0~400	24	6	0~100	
HMS-320	0~400	24	8.5	0~60	
HMS-320T	0~400	24	15	0~60	Roomtemperature+5~50 °C
HMS-330	0~140	22~20	8	0~30	
HMS-340	25~300	0~40	7	0~60	
HMS-350	3000	3	3.7		