

## Water Bath

These baths are available in many different sizes with both digital controls and excellent temperature uniformity. Use optional stainless steel cover for high-temperature applications or cover to speed heating (a Xiangquan Chromtech product).



Model	Heater Power (W)	Volume (L)	Tank Size (mm)	Outside Size (mm)
WB-3100	240	3	235×135×100	265×100×240
WB-5100	350	5	300×150×150	320×175×260
WB-1010	550	10	300×250×150	320×280×260
WB-2010	1100	20	490×295×150	500×300×300